---

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## CONTENTS

---

1	SAFETY PRECAUTIONS	1
2	GENERAL DESCRIPTION	1
2.1	Description	1
2.2	Technical specifications	1
2.3	Warranty	2
3	INSTALLATION	2
3.1	Unpacking the centrifuge	2
3.2	Installation	2
4	OVERVIEW	3
5	OPERATION	4
6	SERVICE AND MAINTENANCE	4
6.1	Centrifuge service	4
6.2	Cleaning the centrifuge	4
6.3	Rotor disinfection	5
7	TECHNICAL SERVICE	6

## 1 SAFETY PRECAUTIONS

---

Before using the *PerfectSpin P* centrifuge for PCR plates for the first time, please read this entire operating manual carefully. To guarantee problem free, safe operation, it is essential to observe the following points:

- NEVER** use the centrifuge in any manner not specified in these instructions.
- NEVER** put hands in the rotor area unless the rotor is completely stopped.
- NEVER** move the centrifuge while the rotor is spinning.
- NEVER** use solvents or flammables near this or other electrical equipment.
- NEVER** centrifuge flammable, explosive or corrosive materials.
- NEVER** centrifuge hazardous materials outside of a hood or proper containment facility.
- NEVER** continue to operate a centrifuge that is vibrating excessively.
- ALWAYS** load the rotor symmetrically. Each plate should be counterbalanced by another plate of the same brand, type and volume.
- ALWAYS** use plates that have been sealed in some manner to avoid spillage.
- ALWAYS** be sure the plates are pushed all the way to the bottom of the rotor slot before operating the centrifuge.
- ALWAYS** locate the centrifuge within easy access to an electrical outlet.

## 2 GENERAL DESCRIPTION

---

### 2.1 Description

The *PerfectSpin P* centrifuge for PCR plates is designed for quick spins of samples in PCR plates. Primarily designed for use with 96 or 384 well plates and small volumes, it will accommodate skirted, non skirted and semi skirted styles.

This manual provides important safety information for the *PerfectSpin P* centrifuge for PCR plates (Cat. No. 91-PSPIN-P / 91-PSPINP-US). It should be kept near the centrifuge for quick and easy reference.

## 2.2 Technical specifications

Model	91-PSPIN-P / 91-PSPINP-US
Max. Speed	2,500 rpm
Max. RCF	500 x g
Max. Rotor Capacity	2 PCR plates
Power Supply	120 V / 50~60 Hz, 0.6 A (Fuse 1.0 AT) 230 V / 50~60 Hz, 0.3 A (Fuse 0.5 AT)
Dimensions (W x H x D)	190 x 183 x 210 mm
Weight (with Rotor)	3 kg

*Note: Specifications are subject to change without prior notice.*

## 2.3 Warranty

PEQLAB guarantees that the *PerfectSpin P* centrifuge for PCR plates you have received has been thoroughly tested and meets its published specification.

However, immediately upon arrival, please check carefully that the shipment is complete and has not been damaged in transit. For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SERVICE'). Please retain all packaging materials until the delivery has been completely checked since this will speed up the return of goods if required and reduce environmental impact. Any form of returns, replacements or credit notes must be agreed in advance by PEQLAB.

This centrifuge has been subjected to thorough testing and quality control. In the unlikely event of a manufacturing fault, our two year warranty (from the date of delivery) covers the centrifuge and the rotor. This warranty becomes invalid in the case of incorrect operation, use of non-standard spare parts or accessories and unauthorized modification of the rotor or centrifuge.

PEQLAB reserves the right to make technical modifications.

## 3 INSTALLATION

### 3.1 Unpacking the centrifuge

Remove the centrifuge from the carton in which it is delivered. All packaging should be retained until it has been established that the centrifuge was proved for transport damage and is working properly.

Confirm that all items are included in the *PerfectSpin P* package:

- *PerfectSpin P* unit with a rotor for two PCR plates
- Instruction manual
- Power cord

For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SERVICE').

### 3.2 Installation

The *PerfectSpin P* centrifuge for PCR plates should be placed on a clean, flat surface with sufficient space around it for ventilation. Be sure that both the bench and the unit's rubber feet are clean. Plug the power cord into the back of the unit and the wall socket. Turn the power switch (on the back of the unit) to the "ON" position.

Mains power and frequency "**must**" be compatible with the information given on the product label on the rear of the centrifuge.

## 4 OVERVIEW

---

This section presents an overview of the *PerfectSpin P* centrifuge's components (see Figure 1) and proper loading of PCR plates (see Figure 2).

Figure 1: *PerfectSpin P* centrifuge for two PCR plates

Figure 2: Proper loading of a PCR plate

## 5 OPERATION

---

To avoid spillage of samples, only plates that have been heat or tape sealed should be placed in the centrifuge. Plates are loaded in the rotor in a vertical position. Lift the clear centrifuge cover to gain access to the rotor (press the "OPEN" button to release the lid if it is locked down). Using the knob in the top of the unit, turn the rotor until the rotor slot lines up with the opening. Carefully slide the first plate into position so that the sealed side faces the center of the rotor and the bottom of the plate faces out. Be sure that the plate is pushed all the way down into the slot. Line up the second slot like the first, using the knob on the top of the unit.

Note: If the rotor is hard to turn, close the lid to start a run to release the friction brake. Load the second plate in the same manner as the first; making sure that it is all the way down in the rotor slot.

**If spinning only one plate, it is extremely important to load a dummy plate of the same brand, style and sample volumes (water can be used) as a balance.**

To begin a run, close the lid of the centrifuge, and the run is started automatically. The centrifuge will automatically accelerate to maximum speed of 2500 rpm, the centrifugal force cannot be adjusted. A run of approximately 15 seconds is sufficient to spin samples to the bottoms of the wells. Pressing the "OPEN" button will engage the friction break and brings the rotor to a stop. Once the rotor has stopped, you can lift the lid completely. To access the plates, turn the knob on the top of the unit until the rotor slot lines up with the opening. If the rotor is hard to turn, close the lid of the centrifuge to start a short run to disengage the friction brake.

Note: A slight vibration may be felt when the centrifuge reaches full speed. Excessive vibration indicates that the rotor is not balanced. DO NOT run the centrifuge unbalanced. Always use a balance plate of the same brand, type and sample volume.

**Attention: Never attempt to operate the centrifuge with a rotor that shows signs of corrosion or mechanical damage. Never centrifuge strongly corrosive materials that may damage the rotors, accessories or bowl of the unit.**

## 6 SERVICE AND MAINTENANCE

---

Always unplug the power plug before service and maintenance.

### 6.1 Centrifuge service

The brushless motor in the *PerfectSpin P* requires no routine maintenance. Any required service should be performed by authorized, qualified personnel only. Repairs performed by unauthorized personnel voids the warranty.

### 6.2 Cleaning the centrifuge

Maintenance of the centrifuge is confined to keeping the rotor and plate channels clean. The outside housing of the centrifuge may be wiped down with a clean, soft cloth. The rotor slots/plate channels may be cleaned with a sponge. Do not pour liquids in or over the centrifuge. Do not immerse the centrifuge in liquid.

Only pH-neutral Detergents with a pH-value of 6-8 may be used for cleaning. Alkaline cleaning agents (pH-value > 8) must not be used! After cleaning, please ensure all parts are dried thoroughly, either by hand or in a hot-air cabinet (maximum temperature +50 °C).

It may also be desired to remove the rotor for cleaning or in the event of spillage. If so continue with the below instruction for disassembly.

### Disassembly:

1. Lay down a soft cloth and turn the unit upside down.
2. Remove the four outer screws.
3. Turn the unit right side up and loosen the rotor knob set screw with a 1/16" allen key and pull of knob.
4. Remove housing by lifting upward.
5. Loosen the rotor connection set screw with a .050" allen key.  
PLEASE NOTE: This rotor contains two set screws directly across from one another. Only the smaller (.050") set screw is to be loosened.
6. The rotor can now be removed for cleaning by lifting upward.

**IMPORTANT:** During reattachment of the rotor, the rotor must be tightened while lifted ½ inch from the baseplate.  
(Typically accomplished by using a ½ inch spacer.)

## 6.3 Rotor disinfection

In case of infectious material spilling into the centrifuge, the rotor and rotor chamber have to be disinfected right after the run. The rotor should be cleaned with a universal, neutral disinfection agent, e.g. in formalin base. A disinfection spray is most suitable in order to easily reach all difficult to access spots.

## 7 TECHNICAL SERVICE

---

Should you have any questions about the *PerfectSpin P*, please visit [www.peqlab.com](http://www.peqlab.com) to find the respective contact person.

Please have the unit's serial number (located on the bottom of the instrument) available when calling. Should an item require return to PEQLAB for service, a decontamination form must be completed first by the user. Items sent without decontamination form will not be accepted.

All returns must be pre-approved by PEQLAB!

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